

Cedex Bio HT analyzer - MAINTENANCE SHEET

Serial No:

Month:

Year:



			Day																																						
	Approximate system time (min)	Approximate operator time (min)	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31								
DAILY																																									
Deproteinize Probes *	6	-																																							
Prime Fluid System *	2	-																																							
WEEKLY																																									
Clean probes and splash guard	2	3																																							
Backup database *	5	-																																							
OPERATORS INITIALS																																									

* automated during BOD - Begin of day